

# Fluorescamine Assay to Determine Protein Content in Polymer Microparticles

From Ainslie Lab @ UNC [ainslielab.web.unc.edu](http://ainslielab.web.unc.edu)

## Reagents

Fluorescamine

Acetone

96 Well Plate (solvent resistant)

## Procedure

- 1) Have samples in 150 ul buffer
- 2) In the first two wells of row A, place 300 ul of the standard
- 3) In the first two wells of the rest of the rows, add 150 ul of sodium acetate buffer
- 4) Do a 2 fold serial dilution down the column, excluding the last row
- 5) Add each sample in 150 ul quantity to two wells (duplicates)
- 6) Determine how much fluorescamine is needed by multiplying the number of wells x 50ul and then add approx. 300 ul extra
- 7) Make the fluorescamine solution in acetone in a 3 mg/ml concentration
- 8) Add 50 ul of fluorescamine solution to each well
- 9) Place film sticker over the plate
- 10) Take plate to plate reader and make sure it is on fluorescence
- 11) Go to excitation and emission. Excitation set at 400nm and Emission at 460nm
- 12) Go to sensitivity, reading should be 40 or above
- 13) Define the wells to read
- 14) Take plate and peel off film, then place in the drawer of reader
- 15) Click on Read, the drawer will close
- 16) Plate will read
- 17) Export file as needed.
- 18) Shut off machine, making sure drawer is closed
- 19) E-mail or save the file and come back to lab
- 20) Open file in Excel
- 21) Make a calibration curve of the standards with the average of the duplicates as the x-axis and the concentration as the y-axis
- 22) Obtain regression equation and analyze results