

<b>Tissue-Culture Vessels</b>	<b>Volume of Media Volume (per item)</b>	<b>Volume of Trypsin (per item)</b>	<b>Effective Growth Area (per item-cm<sup>2</sup>)</b>
<b>Dishes (mm)</b>			
35	2.5 - 3.0 mL	0.2 - 0.3 mL	9.6
60	6.0 - 7.0 mL	0.5 - 0.6 mL	21.3
100	16.0 - 17.5 mL	1.0 mL	58.1
150	45.0 - 50.0 mL	1.5 mL	151.9
<b>Multi-well Plates (wells)</b>			
6	2.5 - 3.0 mL	200 - 300 uL	9.6
12	1.5 - 2.2 mL	100 - 200 uL	3.8
24	0.8 - 1.0 mL	80 - 100 uL	2.0
48	500 - 800 uL	50 - 80 uL	0.75
96	100 - 200 uL	10 - 20 uL	0.32
<b>Flasks</b>			
T-12.5	4.0 - 5.0 mL	250 - 400 uL	12.5
T-25	8.0 - 9.0 mL	500 - 800 uL	25
T-75	20 - 30 mL	1.0 mL	75
T-150	40 - 50 mL	2.0 mL	150
T-175	45 - 55 mL	2.0 mL	175
T-300	150 - 400 mL	4.0 mL	300

Adapted from BD Biosciences Tech Help